

# JOHNSON

## AM-805 Foot Monitor



Three pancake GM detectors configured in a sturdy angled foot rest to make it easy to scan shoes and boots for contamination. Includes a mounted spindle on which to place plastic wrap to cover the detector area and prevent dirt and contamination from falling onto the detectors. Works with any Johnson meter, but especially designed for the RML-200 Lab Counter-top Alarming Ratemeter.

### Specifications

Compatibility	All Johnson survey meters & area monitors
Operating Voltage	900 vdc
Window Density	1.6mg/cm <sup>2</sup>
Range	0-500,000cpm
Dead Time	20 Sec
Energy Response	Alpha--3.5MeV min, Beta--100keV min, Gamma--20keV-2MeV
Sensitive Area	8.5" (21.6cm) X 2" (5cm)
Sensitivity to Cs137	2,500 CPM/mR/hr
Dimensions	9.5" (24cm)L X 11" (28 cm)H X 7" (7.5cm)W
Weight	3.6lbs
Environmental	-28°C--60°C 5%-95% RH non-condensing--IPX4